

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

1 Sample Description

According to this method are examined samples of different types of meat products. Using histochemical staining can demonstrate the presence of starch granules and on their basis to suspect the using of starch or flour, in some cases, illegal production as raw materials in meat products.

2 Detection Limit

The detection limit is 0.1% addition of starch.

3 Time Consumption

Sample preparation
24 hours.

Sample treatment
60 minuts

4 Sampling

4.1 Sample Amount

For the detection it is necessary to take samples in a sufficient amount.

Picked up:

meat product

whole meat product (one package)

4.2 Sampling

Meat products are taken from retail market.

4.3 Sample Fixation

Samples are fixed in 10% water solution of neutral formaldehyde for 24 hours as soon as possible after sampling.

Principles of Fixation:

- samples must be inserted in a special box or in a gauze in a fixative solution as soon as possible,
- sample size should not exceed 1 cm³,

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

- the amount of fixative solution must be at least 20 – 50 times greater than the volume of the fixed sample,
- the sample must be accessible to the fixative solution from all sides, if necessary, it is possible to put the sample on cotton wool or filter paper.

5 Material and Equipment

5.1 Chemicals and Solutions

- alcohol
- acetone
- formol
- hard paraffin grated
- xylene pure
- xylene p. a.
- methylsalicylate
- aluminum sulphate
- nuclear fast red
- indigocarmine
- acid picric
- potassium iodide
- iodine

5.2 Equipment

(uses strictly in accordance of the appropriate manual)

- autotechnikon
- microtom
- fume hood
- staining cuvettes/staining automat
- distillation apparatus
- scales
- pH meter
- thermostat
- hotplates

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

5.3 Laboratory Tools

- gauze
- tweezers
- scalpel
- knife
- cutting mat
- embedding cell
- protective gloves
- markers
- thermometers
- cotton wool
- cotton swabs
- filtration paper
- magnetic stirrer

5.4 Laboratory Glass

- pipettes
- beaker 200 ml
- slides for histochemistry
- cover slips
- funnels

6 Sample Treatment and Preparation

6.1 Preparation of Paraffin Blocks

Samples (5 g) are fixed in 10% water solution of neutral formalin for 24 hours. After fixation, the samples are dehydrated in ascending sequence of alcohol and embedded into paraffin blocks (table 1 and 2).

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

Table 1: Manual preparation of paraffin blocks

| Step No. | Chemicals | Duration |
|----------|-------------------------------|---------------------|
| 1 | Water | 30 min |
| 2 | Alcohol 20% | 30 min. |
| 3 | Alcohol 40% | 60 min. |
| 4 | Alcohol 50% | 60 min. |
| 5 | Alcohol 70% | 120 min. |
| 6 | Alcohol 80% | Per night |
| 7 | Alcohol 96% | 60 min. |
| 8 | Alcohol 100% | 120 min. |
| 9 | Methyl salicylate I | 60 min. |
| 10 | Methyl salicylate II | 60 min. |
| 11 | Methyl salicylate III | Per night |
| 12 | Xylene I | 15 min. |
| 13 | Xylene II | 15 min. |
| 14 | Xylene III | 15 min. |
| 15 | Paraffin I | 120 min. |
| 16 | Paraffin II (p. a.) | 120 min. |
| 17 | Paraffin III (p. a.) | Per night |
| 18 | Embedded into paraffin blocks | 4 th day |

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

Table 2: Preparation of paraffin blocks in autotechnicon

| Step No. | Chemicals | Duration |
|----------|----------------------------|----------|
| 1 | Water | 30 min |
| 2 | Alcohol 50% | 40 min. |
| 3 | Alcohol 70% | 20 min. |
| 4 | Alcohol 96% | 60 min. |
| 5 | Alcohol 96% | 100 min. |
| 6 | Alcohol 96% | 60 min. |
| 7 | Alcohol 100% | 60 min. |
| 8 | Alcohol 100% | 60 min. |
| 9 | Acetone | 20 min. |
| 10 | Xylene I | 20 min. |
| 11 | Xylene II | 20 min. |
| 12 | Paraffin I | 180 min. |
| 13 | Paraffin II (p. a.) | 12 hours |
| 14 | Embed into paraffin blocks | |

6.2 Embedding into Paraffin Blocks

We recommend preparing at least four blocks for each sample. For the purposes of embedding samples into paraffin blocks, commercial medium based on paraffin in combination with bee wax is used. This medium is insoluble in water. Samples saturated with paraffin are embedded in embedding cells by an embedding line. After cooling down, these blocks are prepared for cutting.

6.3 Cutting the Blocks

Tissues embedded in paraffin blocks are cut to 4 μm sections by a microtome. The samples are cut for examination according to the procedure described in the following scheme. The sections are then spread on the water surface and mounted to slides.

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

Cutting Scheme

block A

section 1 – cut off 50 µm – **section 2** – cut off 50 µm – **section 3** – cut off 50 µm – **section 4** – cut off 50 µm – **section 5** – cut off 50 µm – **section 6**

block B

section 7 – cut off 50 µm – **section 8** – cut off 50 µm – **section 9** – cut off 50 µm – **section 10** – cut off 50 µm – **section 11** – cut off 50 µm – **section 12**

block C

section 13 – cut off 50 µm – **section 14** – cut off 50 µm – **section 15** – cut off 50 µm – **section 16** – cut off 50 µm – **section 17** – cut off 50 µm – **section 18**

block D

section 19 – cut off 50 µm – **section 20** – cut off 50 µm – **section 21** – cut off 50 µm – **section 22** – cut off 50 µm – **section 23** – cut off 50 µm – **section 24**

6.4 Stainig

Sectiones are necessary before staining get rid of paraffin. Paraffin sections dissolves in the solvent (xylene) and through ethanol are the incisions converted into water because histological dyes are mostly soluble in water. We can dyed by hand in special cuvettes or in the staining automat. Dyes and dye mixtures are applied according to the procedure described below.

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

Lugol Calleja

Objective: highlight the starch granules

Procedure:

| | | |
|-------------------------|--------------------------------------|---------|
| Getting rid of paraffin | xylene..... | 10 min. |
| | 100 % alcohol + ether (2/3+1/3)..... | 10 min. |
| Nuclear fast red..... | | 15 min. |
| Distilled water..... | | Wash |
| Lugol's solution..... | | 5 min. |
| Distilled water..... | | Wash |
| Solution B Calleja..... | | 5 min. |
| Distilled water..... | | Wash |
| Dewatering (alcohol) | 96%..... | Wash |
| | 100%..... | Wash |
| xylene I (pure) | | 5 min. |
| xylene II (p. a.) | | 5 min. |

Preparation of Solutions:

Nuclear fast red:

| | |
|--------|-------------------|
| 10 g | aluminum sulphate |
| 0,1 g | nuclear fast red |
| 100 ml | distilled water |

- dissolve over low heat, simmer briefly. After chilling filtered.

Solution B Calleja:

| | |
|--------|---|
| 100 ml | indigocarmine 1% (1,0 g indigocarmine+ 100 ml distilled water) |
| 200 ml | acid picric (saturated solution) |

Lugol's solution:

| | |
|--------|------------------|
| 300 ml | distilled water |
| 2,0 g | potassium iodide |
| 1,0 g | iodide |

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

Stainig results:

| | |
|--|----------|
| core | – red |
| muscle | – green |
| elastic connective tissue | – yellow |
| fibrous connective tissue, bone, cartilage | – blue |
| starch grain | – brown |

6.5 Mounting of Stained Slides

Stained sections are mounted between a slide and a cover slip with a suitable mounting medium which does not interfere with the color. Usually, used synthetic resins are insoluble in water, so it is necessary to dewater the sections in an ascending sequence of alcohol and xylene.

Process for Manual Mounting:

1. put a drop of mounting medium on the margin of the slide.
2. place a cover slip on the margin of the slide under the angle of 45°.
3. carefully and slowly move the cover glass in order to avoid formation of bubbles.
4. put slides in thermostat with 60°C for the night.
5. carefully clean the edges using alcohol and a razor.

It is also possible to use a mounting automate.

7 Microscopic Examination and Evaluation of Results

The stained sections are examined by the light microscope with a lower magnification (e.g. 32x or 40x), for the study of detail is used higher magnification. Usually, only qualitative examination is performed. Describes the presence of different types of tissues in the examined samples. It is possible to focus only on identification of selected tissues highlighted by special staining. The identification of tissues of animal and vegetable origin must be based on data from the literature. For comparison use samples prepared in the laboratory and also the schematic pictures and photos from the literature.

8 Documentation

It is recommended to make a laboratory protocol for each sample with the following information:

- number of sample
- date of sampling
- type of product
- producer

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

- ingredients (if known)
- sample treatment and preparation (fixation, staining, etc.)
- examination results


For clear identification, it is also recommended to label sample container by sample name and number (the same as written in the protocol). Blocks and slides should be labeled in the same way. It is possible to label samples by bar-codes.

All samples with photo documentation are archived.

9 Results

Starch granules are stained brown. Due to the loss characteristics can not determine the type of starch grains.

10 Photo Documentation



Starch granules in a meat product – “Kabanos”, staining Lugol - Calleja, magnification 200x

| Protocol 05 | |
|-----------------------------|--|
| Monitored ingredient | Starch grains |
| Foodstuff | Meat products |
| Examination | Histochemistry, Lugol Calleja stainig |
| Short protocol/full version | Full version |

11 List of Abbreviations

LC – Lugol – Calleja staining